

RadiantClariX - AI X-Ray Analysis Application
User Manual & Documentation

Version: 1.0.0
Release Date: October 2025
Last Updated: October 30, 2025

Table of Contents

1. Introduction
2. Key Features
3. How to Start the Application
4. User Guide
5. Settings & Features
6. Troubleshooting
7. Technical Specifications
8. Privacy & Security
9. Support
10. Important Disclaimers
11. About the AI Models
12. Version Information

Introduction

RadiantClariX is an advanced AI-powered medical imaging application that analyzes X-ray images and provides detailed diagnostic reports. The app uses machine learning models to detect abnormalities in chest X-rays and identify bone fractures with high accuracy.

Perfect for:

- Healthcare professionals seeking quick preliminary analysis
- Medical students learning radiology
- Clinics needing AI-assisted diagnostics
- Research and educational purposes

Key Features

Dual AI Models:

1. Chest X-Ray Analyzer - Generates detailed textual descriptions of chest X-ray images. Describes visible anatomical features, abnormalities, and medical devices. Provides both technical medical terminology and plain language explanations. Uses BLIP AI technology for image captioning

2. Bone Fracture Detector - Identifies seven bone injury types with bounding boxes and confidence scores. Detects 7 types of bone injuries:

User-Friendly Interface:

- Dark & Light Theme
- Easy Image Upload
- Real-time Analysis
- Scan History
- PDF Export

Secure & Private:

- Authentication system
- Encrypted data storage
- Option to delete accounts and data

How to Start the Application

Step 1: Start the Backend Services

1. Open Project file in VSCode or in any other IDE
2. Navigate to the backend folder:
`cd D:\xray-app\radiantclarix\backend`
3. Type “**.\start-all-services.ps1**” and hit ENTER.

Step 2: Start the Frontend App

1. Open Project file in VSCode or in any other IDE
2. Navigate to app folder:
`cd D:\xray-app\radiantclarix`
3. Type: “**npm start**” and hit ENTER

Step 3: Open the App

Use Expo Go on mobile or emulator/web browser.

User Guide

Registration & Login:

- Register or Login using email and password.

Dashboard:

- Upload & Analyze, View History, Settings, Logout.

Upload & Analyze:

- Select model (Chest/Bone)
- Enter patient name
- Upload image (camera or gallery)
- Analysis starts automatically.

Understanding Results:

- Technical Report (for professionals)
- Plain Language Report (for patients)
- Generate PDF, Start New Scan, View Annotated Image.

Scan History:

- Full record of all scans with reports and dates.

Settings Page:

- Account Info, Theme, Change Password, Delete History, Delete Account, Logout.

Troubleshooting

App Won't Start:

- Ensure backend services are running.
- Restart both backend and frontend.

Login Issues:

- Check credentials and internet connection.

Upload Issues:

- Grant permissions, restart app, ensure models running.

History Not Updating:

- Refresh or restart app.

Performance:

- Restart backend, use smaller images, clear cache.

Technical Specifications

System Requirements:

Backend - Windows 10/11, Python 3.8+, Node.js 14+, 8GB RAM, 5GB disk.

Frontend - Android/iOS device with Expo Go or browser.

Ports Used:

5000 - Backend

8502 - Chest X-Ray AI Model

8503 - Bone Fracture Model

8081 - Expo Dev Server

Privacy & Security

Data Protection:

- Passwords encrypted and hashed
- MongoDB secure storage

User Rights:

- View, delete, or export data anytime.

Recommendations:

- Use strong passwords
- Logout on shared devices
- Backup reports regularly

Support

Q: Is this suitable for medical diagnosis?

A: No, use for preliminary analysis only.

Q: Offline mode?

A: Internet required.

Q: How many scans can I store?

A: Unlimited. Reports will be available for 15 days only from the time they are generated.

Q: Can I share reports?

A: Yes, via PDF export.

Q: Data privacy?

A: Fully private and user-only access.

Important Disclaimers

Medical Disclaimer:

- Educational use only.
- Do not replace professional diagnosis.
- Consult professionals for treatment.

Accuracy Notice:

- AI confidence scores vary.
- False positives/negatives possible.
- Models updated regularly for improved accuracy.

About the AI Models

Chest X-Ray Model:

- BLIP (Bootstrapping Language-Image Pretraining)
- Detects lung/heart issues, medical devices.

Bone Fracture Model:

- Faster R-CNN with ResNet-50 backbone
- Detects seven fracture types with bounding boxes and confidence scores.

Version Information

Application Version: 1.0.0

Release Date: October 2025

Last Updated: October 30, 2025

Thank you for using RadiantClariX.
For technical support, contact your system administrator.